

TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

24.0cm (9.6")



SKU: 051-080-24

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM078
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-24
AUTHORS	HEGAR , MAYO
Medicon	10.18.74
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM078

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-24
AUTHORS	HEGAR, MAYO
Medicon	10.18.74
SPECIFICATIONS	serrated , straight