


TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

24.0cm (9.6")



SKU: 051-080-24

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM078 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-642-24 |
| AUTHORS | HEGAR , MAYO |
| Medicon | 10.18.74 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| Aesculap | BM078 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-642-24 |
| AUTHORS | HEGAR, MAYO |
| Medicon | 10.18.74 |
| SPECIFICATIONS | serrated , straight |