

TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-080-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-26
AUTHORS	HEGAR , MAYO
Medicon	10.18.76
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-26

AUTHORS

[HEGAR](#), [MAYO](#)

Medicon

10.18.76

SPECIFICATIONS

[serrated](#), [straight](#)