


TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

30.0cm (12")



SKU: 051-080-30

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM069
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-30
AUTHORS	HEGAR , MAYO
Medicon	10.18.80
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM069

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-30
AUTHORS	HEGAR, MAYO
Medicon	10.18.80
SPECIFICATIONS	serrated , straight