


## TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

30.0cm (12")



**SKU:** 051-080-30

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM069
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-642-30
<b>AUTHORS</b>	<a href="#">HEGAR</a> , <a href="#">MAYO</a>
<b>Medicon</b>	10.18.80
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM069

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-642-30
<b>AUTHORS</b>	<a href="#">HEGAR, MAYO</a>
<b>Medicon</b>	10.18.80
<b>SPECIFICATIONS</b>	<a href="#">serrated, straight</a>