

TC MAYO-HEGAR NEEDLE HOLDER LH

straight; serrated; for left hand use

18.0cm (7")



SKU: 051-080-18/LH

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR , MAYO
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	HEGAR , MAYO
SPECIFICATIONS	serrated , straight