

TC MAYO-HEGAR NEEDLE HOLDER LH

straight; serrated; for left hand use

18.0cm (7")




SKU: 051-080-18/LH

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | HEGAR , MAYO |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | HEGAR , MAYO |
| SPECIFICATIONS | serrated , straight |