

TC MAYO-HEGAR NEEDLE HOLDER

straight; serrated

16.0cm (6.4")

SKU: 051-080-16

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM065
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-16
AUTHORS	HEGAR , MAYO
Medicon	10.18.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM065

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-642-16
AUTHORS	HEGAR, MAYO
Medicon	10.18.65
SPECIFICATIONS	serrated, straight