

TC MICRO-ADSON DRESSING FORCEPS

straight; serrated; jaw 1,1+/-0.2

15.0cm (6")

SKU: 021-076-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CUSHING
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	CUSHING
SPECIFICATIONS	serrated , straight