

TC MICRO-ADSON DRESSING FORCEPS

straight; serrated; jaw 1,1+/-0.2

12.0cm (4.8")

SKU: 021-076-12

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	ADSON , MICRO
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
AUTHORS	ADSON , MICRO
SPECIFICATIONS	serrated , straight