


TC MICRO RYDER NEEDLE HOLDER

straight; serrated

15.0cm (6")



SKU: 051-132-15

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM055
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-15
AUTHORS	MICRO , RYDER
Medicon	10.28.65
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM055

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-15
AUTHORS	MICRO, RYDER
Medicon	10.28.65
SPECIFICATIONS	serrated, straight