


TC MICRO RYDER NEEDLE HOLDER

straight; serrated

18.0cm (7.2")

SKU: 051-132-18

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM056
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-18
AUTHORS	MICRO , RYDER
Medicon	10.28.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM056

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-18
AUTHORS	MICRO, RYDER
Medicon	10.28.68
SPECIFICATIONS	serrated, straight