


TC MICRO RYDER NEEDLE HOLDER

straight; serrated

23.0cm (9.2")



SKU: 051-132-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-23
AUTHORS	MICRO , RYDER
Medicon	10.28.73
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-23

AUTHORS

[MICRO](#), [RYDER](#)

Medicon

10.28.73

SPECIFICATIONS

[serrated](#), [straight](#)