


TC MICRO RYDER NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-132-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-26
AUTHORS	MICRO , RYDER
Medicon	10.27.75
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-654-26

AUTHORS

[MICRO](#), [RYDER](#)

Medicon

10.27.75

SPECIFICATIONS

[serrated](#), [straight](#)