


## TC MICRO RYDER NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



**SKU:** 051-132-13

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM054
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-654-13
<b>AUTHORS</b>	<a href="#">MICRO</a> , <a href="#">RYDER</a>
<b>Medicon</b>	10.28.63
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM054

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-654-13
<b>AUTHORS</b>	<a href="#">MICRO, RYDER</a>
<b>Medicon</b>	10.28.63
<b>SPECIFICATIONS</b>	<a href="#">serrated, straight</a>