


TC MICROVASCULAR NEEDLE HOLDER

straight; serrated

23.0cm (9.2")



SKU: 051-122-23

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-666-04 |
| AUTHORS | MICROVASCULAR |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 20-666-04 |
| AUTHORS | MICROVASCULAR |

SPECIFICATIONS

[serrated](#), [straight](#)