

## TC NEEDLE HELDER; STRAIGHT

shank handle with ratchet

Luer Lock; Ø 5.0mm; 450mm



**SKU:** 614-005-54

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)