

TC NEEDLE HOLDER Ø 3MM, CURVED LEFT

with Aluminium handle, teflon coated

WL 290mm, LL-Connection



SKU: 614-110-33

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#), [Pediatric Laparoscopy](#)

TYPE [LOFTLINE](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#), [Pediatric Laparoscopy](#)

TYPE [LOFTLINE](#)