

## TC NEEDLE HOLDER Ø 3MM, STRAIGHT

with Aluminium handle, teflon coated

WL 290mm, LL-Connection



**SKU:** 614-100-33

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#), [Pediatric Laparoscopy](#)

**TYPE** [LOFTLINE](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#), [Pediatric Laparoscopy](#)

**TYPE** [LOFTLINE](#)