

## TC NEEDLE HOLDER Ø 5MM, CONVEX/ CONCAVE

with Aluminium-pistol handle

teflon coated, WL 330mm, LL-Connection



**SKU:** 614-520-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [TROMPLINE](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [TROMPLINE](#)