

TC NEEDLE HOLDER Ø 5MM, CURVED RIGHT

with Aluminium handle, teflon coated

WL 450mm, LL-Connection



SKU: 614-115-54

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Needle Holders
TYPE	LOFTLINE

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Needle Holders
TYPE	LOFTLINE