

TC NEEDLE HOLDER Ø 5MM, "S-SHAPED"

with Aluminium-pistol handle

teflon coated, WL 330mm, LL-Connection



SKU: 614-535-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [TROMPLINE](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [TROMPLINE](#)