

TC OEHLER TISSUE FORCEPS

straight; serrated; 1X2 teeth

30.0cm (12")

SKU: 021-228-30

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-260-30
AUTHORS	SEMKEN
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-260-30
AUTHORS	SEMKEN

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)