


TC OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

19.0cm (7.6")




SKU: 051-310-19

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-672-18
AUTHORS	HEGAR , OLSEN
Medicon	10.60.69
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-672-18

AUTHORS

[HEGAR](#), [OLSEN](#)

Medicon

10.60.69

SPECIFICATIONS

[serrated](#), [straight](#)