


TC OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

14.0cm (5.6")



SKU: 051-310-14

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM128
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-672-14
AUTHORS	HEGAR , OLSEN
Medicon	10.60.64
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM128

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-672-14
AUTHORS	HEGAR, OLSEN
Medicon	10.60.64
SPECIFICATIONS	serrated , straight