

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

20.0cm (8")


SKU: 021-026-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD157 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 12-258-21 |
| AUTHORS | POTTS , SMITH |
| Medicon | 06.02.70 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD157 |

| | |
|-----------------------|---|
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 12-258-21 |
| AUTHORS | POTTS, SMITH |
| Medicon | 06.02.70 |
| SPECIFICATIONS | serrated , straight |