

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

20.0cm (8")

SKU: 021-026-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD157
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-21
AUTHORS	POTTS , SMITH
Medicon	06.02.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD157

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-21
AUTHORS	POTTS, SMITH
Medicon	06.02.70
SPECIFICATIONS	serrated , straight