

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

23.0cm (9.2")

SKU: 021-026-23

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-23
AUTHORS	POTTS , SMITH
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-23
AUTHORS	POTTS , SMITH

SPECIFICATIONS

[serrated](#), [straight](#)