

## TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

23.0cm (9.2")

**SKU:** 021-026-23

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-258-23
<b>AUTHORS</b>	<a href="#">POTTS</a> , <a href="#">SMITH</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-258-23
<b>AUTHORS</b>	<a href="#">POTTS</a> , <a href="#">SMITH</a>

## SPECIFICATIONS

[serrated](#), [straight](#)