

TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

25.0cm (10")

SKU: 021-026-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD159
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-258-25
AUTHORS	SEMKEN
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD159
MATERIAL	Tungsten Carbide(TC)

KLS Martin

12-258-25

AUTHORS

[SEMKEN](#)

SPECIFICATIONS

[serrated](#), [straight](#)