

## TC POTTS-SMITH DRESSING FORCEPS

straight; serrated

25.0cm (10")


**SKU:** 021-026-25

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                              |
| <b>TYPE</b>           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b>       | BD159                                                             |
| <b>MATERIAL</b>       | <a href="#">Tungsten Carbide(TC)</a>                              |
| <b>KLS Martin</b>     | 12-258-25                                                         |
| <b>AUTHORS</b>        | <a href="#">SEMKEN</a>                                            |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a>               |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>                                              |
| <b>TYPE</b>     | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b> | BD159                                                             |
| <b>MATERIAL</b> | <a href="#">Tungsten Carbide(TC)</a>                              |

**KLS Martin**

12-258-25

**AUTHORS**

[SEMKEN](#)

**SPECIFICATIONS**

[serrated](#), [straight](#)