


TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

18.0cm (7.2")



SKU: 051-134-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM045
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-18
AUTHORS	RYDER
Medicon	10.26.68
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM045

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-18
AUTHORS	RYDER
Medicon	10.26.68
SPECIFICATIONS	serrated , straight