


TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-134-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM046
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-20
AUTHORS	RYDER
Medicon	10.26.70
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM046

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-20
AUTHORS	RYDER
Medicon	10.26.70
SPECIFICATIONS	serrated , straight