


## TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

23.0cm (9.2")



**SKU:** 051-134-23

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-650-23
<b>AUTHORS</b>	<a href="#">RYDER</a>
<b>Medicon</b>	10.26.73
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	20-650-23

## AUTHORS

[RYDER](#)

## Medicon

10.26.73

## SPECIFICATIONS

[serrated](#), [straight](#)