


TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

26.0cm (10.4")



SKU: 051-134-26

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM047
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-26
AUTHORS	RYDER
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM047
MATERIAL	Tungsten Carbide(TC)

KLS Martin

20-650-26

AUTHORS

[RYDER](#)

SPECIFICATIONS

[serrated](#), [straight](#)