


TC RYDER "FRENCH EYE" NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



SKU: 051-134-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM044
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-13
AUTHORS	RYDER
Medicon	10.26.63
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM044

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-650-13
AUTHORS	RYDER
Medicon	10.26.63
SPECIFICATIONS	serrated , straight