


TC SAROT NEEDLE HOLDER

straight; serrated

20.0cm (8")




SKU: 051-224-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	SAROT
Medicon	10.22.72
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	SAROT
Medicon	10.22.72

SPECIFICATIONS

[serrated](#), [straight](#)