


## TC SAROT NEEDLE HOLDER

straight; serrated

26.0cm (10.4")




**SKU:** 051-224-26

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM023
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	26-640-26
<b>AUTHORS</b>	<a href="#">SAROT</a>
<b>Medicon</b>	10.22.77
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM023

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	26-640-26
<b>AUTHORS</b>	<a href="#">SAROT</a>
<b>Medicon</b>	10.22.77
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>