

TC SC CC HEGEMANN-DIETRICH VASCULAR SCI

45° curved

19.0cm (7.6")




SKU: 014-550-45

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Scissors for Vascular Surgery |
| MATERIAL | Ceramic(CC) , Super-Cut(SC) , Tungsten Carbide(TC) |
| AUTHORS | DIETRICH , HEGEMANN |
| SPECIFICATIONS | curved |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Scissors for Vascular Surgery |
| MATERIAL | Ceramic(CC) , Super-Cut(SC) , Tungsten Carbide(TC) |
| AUTHORS | DIETRICH , HEGEMANN |
| SPECIFICATIONS | curved |