

TC SEMKEN DRESSING FORCEPS

straight; serrated

15.0cm (6")


SKU: 021-040-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-264-15
AUTHORS	SEMKEN
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-264-15
AUTHORS	SEMKEN

SPECIFICATIONS

[serrated](#), [straight](#)