

## TC SEMKEN TISSUE FORCEPS

straight; serrated; 1X2 teeth

15.0cm (6")


**SKU:** 021-240-15

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-240-16
<b>AUTHORS</b>	<a href="#">SEMKEN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-240-16
<b>AUTHORS</b>	<a href="#">SEMKEN</a>

## SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)