


TC SWEDISH MODELL NEEDLE HOLDER

straight; serrated

17.0cm (6.8")



SKU: 051-136-17

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-652-17
AUTHORS	SCHWEDISCHES MODELL
Medicon	10.39.67
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-652-17

AUTHORS

[SCHWEDISCHES MODELL](#)

Medicon

10.39.67

SPECIFICATIONS

[serrated](#), [straight](#)