


TC SWEDISH MODELL NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-136-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | SCHWEDISCHES MODELL |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Needle Holders |
| MATERIAL | Tungsten Carbide(TC) |
| AUTHORS | SCHWEDISCHES MODELL |
| SPECIFICATIONS | serrated , straight |