


TC SWEDISH MODELL NEEDLE HOLDER

straight; serrated

20.0cm (8")



SKU: 051-136-20

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	SCHWEDISCHES MODELL
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
MATERIAL	Tungsten Carbide(TC)
AUTHORS	SCHWEDISCHES MODELL
SPECIFICATIONS	serrated , straight