

TC WANGENSTEEN DRESSING FORCEPS

straight; cross serrated

23.0cm (9.2")

SKU: 021-140-23

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD173
MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-266-23
AUTHORS	MC INDOE
Medicon	07.16.73
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD173

MATERIAL	Tungsten Carbide(TC)
KLS Martin	12-266-23
AUTHORS	MC INDOE
Medicon	07.16.73
SPECIFICATIONS	serrated , straight