

## TC WANGENSTEEN DRESSING FORCEPS

straight; cross serrated

15.0cm (6")

**SKU:** 021-140-15

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD171
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-266-15
<b>AUTHORS</b>	<a href="#">WANGENSTEEN</a>
<b>Medicon</b>	07.16.65
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD171

<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	12-266-15
<b>AUTHORS</b>	<a href="#">WANGENSTEEN</a>
<b>Medicon</b>	07.16.65
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>