

TC WEBSTER NADELHALTER

straight; smooth

13.0cm (5.2")



SKU: 051-028-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM014
MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-630-13
AUTHORS	WEBSTER
Medicon	10.10.62
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM014

MATERIAL	Tungsten Carbide(TC)
KLS Martin	20-630-13
AUTHORS	WEBSTER
Medicon	10.10.62
SPECIFICATIONS	straight