

TC WIRE CUTTING PLIER; CUTTING WIDTH 18M

soft wire up to Ø2.0mm; hard to Ø1.2mm

double convertible; 15.0cm (6")



SKU: 190-770-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX150
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-524-15
Medicon	14.30.65

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX150
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-524-15

Medicon

14.30.65