

## TC WIRE CUTTING PLIER

soft wire up to Ø3.0mm; hard to Ø2.2mm

double convertible; 23.0cm (9.2")


**SKU:** 190-774-23

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Forceps</a>
<b>Aesculap</b>	LX156
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-523-24
<b>Medicon</b>	76.45.73



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Forceps</a>
<b>Aesculap</b>	LX156
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-523-24

**Medicon**

76.45.73