

TC WIRE CUTTING PLIER; CUTTING WIDTH 10MM

soft wire up to Ø2.0mm; hard to Ø1.5mm

double convertible; 18.0cm (7.2")



SKU: 190-772-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX158
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-526-18
Medicon	14.30.28

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX158
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-526-18

Medicon

14.30.28